

## PROTOCOL: CHANNELS SEPARATION + OMEFILES (DATA FROM SPINNING DISK)


The following protocol uses the functions `ch_separation02.m` and `tiffs2ome01.m` that are in 'labo5\paola\ch\_separation\functions'. Use this protocol and these functions if:

- YOUR DATA ARE FROM THE SPINNING DISK
- THE CHANNELS ARE NOT SPLITTED
- YOU HAVE ONLY 1 DEVELOPMENTAL TIME (NO TIMELAPSES)

Start Matlab


Copy the folder `Labo5\Paola\ch_separation\functions` on your desktop

Choose the folder you created as Matlab directory, by clicking here:



And choosing the right folder.


A list of functions will appear.



Double click on **ch\_separation02.m**

Define as **project\_path** the path of the folder where the data from the spinning disk are saved (ONLY ONE EXPERIMENT IN ONE FOLDER!).

The path has to be between ' and \;



Define the number of channel as **nC**, the number of frames as **nt** and the number of z planes as **nz**.  
Don't change nr!

<pre> % number of channels nC=2 % number of repeats (developmental time points) nr=1 % number of frames nt=200 % number of z planes nz=201 dz=1  %% define order of the files depending on how they were exported % order as exported from sp8 'zrct' % order='zrct'; % FromLightSheet=1;% set to 1 for 16 bit images  %% file exported from fiji (single developmental point) exported from sp8 'zrct' % order='tzc';  %% file exported from lightSheet microscope (single developmental point) exported from sp8 'zrc' order='tzc'; % only single channel tested so far %% voxelSizeY=1; </pre>	
---	--

Run the code.

A folder named **Separated** will appear in the same folder where the data are saved.

Inside a list of images is saved in the order tzc.

A folder named **Omefiles** will appear in the folder **Separated**.  
Inside a list of ome files is saved in the order rz.  
The ome files are the files you need for the realignment!